ESPAÑOL 15-28

DEUTSCHE 29-42



FRANÇAIS 43-56

ITALIANO 57-70

Pulse Oximeter AE174

QUICK START GUIDE



SCREEN DISPLAY

OPERATION

Install Battery



- ▶ Install two AAA batteries into the battery compartment.

OPERATION

Using the Lanyard

- Tie the lanyard through the hole in the rear of the Pulse Oximeter.

To Operate

- Place one of your fingers into the Pulse Oximeter to the end and press the switch button one time on front panel to turn it on.

Note: Keep your hands still for the reading.

SETTING

V2.12
Alm setup *
Alm on
Beep off
Demo off
Restore
Brightness 4
Exit

Alm setup	To enter the Alarm Setup submenu
Alm	To turn audible alarm on/off
Beep	To turn audible beep for pulse beat on/off
Demo	To demonstrate screen display without real measurement
Restore	To restore default settings
Brightness	To adjust screen brightness from 1-10. Level 10 is the brightest. Default level is 4
Exit	To return to the measuring interface

SETTING

V2.12
Sounds Setup *
SpO₂ Alm Hi 100
SpO₂ Alm Lo 94
PR Alm Hi 130
PR Alm Lo 50
+/- +
Exit

SpO ₂ Alm Hi	Upper limit for SpO ₂ alarm
SpO ₂ Alm Lo	Lower limit for SpO ₂ alarm
PR Alm Hi	Upper limit for PR alarm
PR Alm Lo	Lower limit for PR alarm
+/-	Select + to increase value. Select - to decrease value.
Exit	To return to the measuring interface

PULSE SOUND

- When the Oximeter is ON, press and hold POWER (⏻) button until the setting menu.
- Short press (<1 sec) POWER (⏻) button to scroll down the “*” symbol to “Beep”.
- Long press (>1 sec) POWER (⏻) button to change “Beep” to ON or OFF.

V2.12
Alm setup
Alm
Beep * on
Demo off
Restore off
Brightness 4
Exit

BRIGHTNESS ADJUSTMENT

- Under the setting menu.
- Short press (<1 sec) POWER (⏻) button to scroll down the “*” symbol to “Brightness”.
- Long press (>1 sec) POWER (⏻) button to adjust brightness from 1 (Low) to 5 (High).


A screenshot of a black screen with white and yellow text. The text is arranged in a list-like format. The first column contains the following items: V2.12, Alm setup, Alm, Beep, Demo, Restore, and Brightness (highlighted in yellow). To the right of these items, there are two columns of text: 'on' and 'off' (aligned with 'Alm'), and 'off' (aligned with 'Demo'). At the bottom right, there is a yellow asterisk (*) and the number '4'. At the bottom center, the word 'Exit' is displayed in white.

V2.12		
Alm setup		
Alm	on	
Beep	off	
Demo	off	
Restore		
Brightness		
	Exit	* 4

OUT-OF RANGE ALARM

- When SpO₂ or PR reading exceeds a defined range, an alarm will trigger and the SpO₂ or PR reading will start to flash respectively. In addition, an audible alarm can be turned on by setting the “Alm” under the setting menu to ON.

V2.12
Alm setup
Alm * on
Beep off
Demo off
Restore off
Brightness 4
Exit 4

- Short press (<1 sec) POWER (⏻) button to move the “*” symbol to “Alm” option.
- Long press (>1 sec) POWER (⏻) button to turn alarm function ON or OFF.

OUT-OF RANGE ALARM

- The upper and lower limit for the SpO₂ and PR alarm can be user defined. Under the setting menu, select “Alm Setup” long press (>1 sec) POWER (⏻) button to enter the Alarm Setup submenu.
- To increase a value, ‘+’ should be selected.
- To decrease a value, ‘-’ should be selected.
- Short press (<1 sec) POWER (⏻) button to move the “*” symbol to “SpO₂ Alm Hi”.
- Long press (>1 sec) POWER (⏻) button to change to your desired value. Follow the same steps to change other values.

V2.12
Sounds Setup
SpO₂ Alm Hi * 100
SpO₂ Alm Lo 94
PR Alm Hi 130
PR Alm Lo 50
+/- +

Exit

NOTE

Measurements may be inaccurate in the following situations:

- In a very cold environment or cold fingers
- Excessive finger movement during measurement
- When finger isn't placed properly into the device

Excessive ambient light, such as fluorescent lamp or direct sunlight, may affect the result.

Enamel or acrylic fingernail polish or other fingernail applications may give inaccurate readings.

FAQ

What is Pulse Oximeter?

A pulse oximeter is a non-invasive device that indirectly monitors blood oxygen saturation (SpO_2) and pulse rate (heart rate). Pulse oximeters provide an easy way of assessing your blood oxygen level and pulse rate.

What is SpO_2 ?

SpO_2 is also known as oxygen saturation. Oxygen saturation is a measure of how much oxygen the blood is carrying as a percentage of the maximum it could carry.

What kind of conditions may cause an inaccurate reading?

Cold hands, poor circulation, very weak pulse, movement, fingernail polish and acrylic nails may cause inaccurate results.

FAQ

What is the normal range for pulse rate?

The normal resting range for pulse rate is typically considered from 60-100 beats per minute. Ask your health professional this question as it pertains to you.

What is the normal range of SpO₂?

The normal range for SpO₂ is typically considered from 95% - 99%. The SpO₂ measurement may be lower for people who live at high altitudes. Ask your health professional this question as it pertains to you.

The SpO₂ is not changing-it's stuck?

SpO₂ does not change like pulse rate. It is slow to change.

The pulse rate is changing rapidly.

Your heart rate changes with emotions, excitement and exercise.

EXTENDED WARRANTY

Go to www.AlcedoHealth.com/Warranty to register your product and receive one FREE extra year warranty.

Please read the instruction manual for complete product information.

Have questions?

DO NOT RETURN THIS PRODUCT.

Call our Toll-Free Consumer Help Line at **833-682-8655** or email support.eu@alcedohealth.com for assistance with your Alcedo products.

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

A screenshot of a device's configuration menu. The menu is displayed on a black background with white and yellow text. The items listed are: V2.12, Alm setup, Alm, Beep, Demo, Restore, and Brightness. To the right of 'Beep' is an asterisk (*), and to the right of 'Demo' are the words 'on' and 'off' stacked vertically. At the bottom right, the number '4' is displayed. The word 'Exit' is at the bottom center in yellow.

V2.12	
Alm setup	
Alm	
Beep	*
Demo	on
Restore	off
Brightness	4
Exit	

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
V2.12
Alm setup
Alm
Beep
Demo
Restore
Brightness
Exit

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V2.12		
Alm setup		
Alm	*	on
Beep		off
Demo		off
Restore		
Brightness		4
Exit		

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V2.12		
Alm setup		
Alm		
Beep	*	on
Demo		off
Restore		off
Brightness		4
Exit		

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
A screenshot of a black rectangular menu box with white and yellow text. The menu items are listed vertically on the left: 'V2.12', 'Alm setup', 'Alm', 'Beep', 'Demo', 'Restore', 'Brightness' (in yellow), and 'Exit' (in yellow). To the right of 'Alm' is an asterisk '*' and the words 'on', 'off', 'off' stacked vertically. To the right of 'Brightness' is the number '4'.

V2.12	
Alm setup	
Alm	*
Beep	on
Demo	off
Restore	off
Brightness	4
Exit	

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